

Laboratory Detection and Reporting of Carbapenem-Resistant Enterobacteriaceae (CRE)

**CLSI Outreach Working Group
Spring, 2016**

After review of this program, you will be able to:

- ◆ List current **CLSI recommendations** for antimicrobial susceptibility testing (AST) and reporting of **CRE**.
- ◆ Describe the differences between **CRE** (Carbapenem-resistant Enterobacteriaceae) and **CPE** (Carbapenemase-producing Enterobacteriaceae).
- ◆ Explain the significance of CRE and CPE from **clinical and epidemiological perspectives**.
- ◆ Describe when and where to go for **help** when physicians need more information than your laboratory can provide for a potential CRE isolate.

Acronyms Used in this Presentation

- ◆ **CRE** = Carbapenem-R Enterobacteriaceae
- ◆ **CPE** = Carbapenemase-Producing Enterobacteriaceae
- ◆ Also **CRO, CRKP, CPKP, CPO**
 - “O” = organism
 - “KP” = *Klebsiella pneumoniae*

β -Lactams and GNRs


Basic Concepts

Common β -Lactam Agents Active Against Gram-Negative Rods (GNR)


β -lactam

<u>Penicillins</u>	<u>Inhibitor Combos</u>	<u>Cephalosporins</u>	<u>Cephamycin</u>	<u>Carbapenems</u>
Amoxicillin	Amoxicillin-clav	Cefazolin (1)	Cefoxitin	Doripenem
Ampicillin	Ampicillin- sulb	Cefuroxime (2)	Cefotetan	Ertapenem
Piperacillin	Piperacillin-tazo	Cefotaxime (3)	<u>Monobactam</u>	Imipenem
Ticarcillin	Ticarcillin-clav	Ceftazidime (3)	Aztreonam	Meropenem
	Ceftolozane-tazo	Ceftriaxone (3)		
	Ceftaz-avibactam	Cefepime (4)		

β -Lactams




β -Lactam
(Penicillin)




β -Lactam
(Penicillin)

All β -lactams contain a 4-membered β -lactam ring that is essential for activity


Different β -Lactams




Penicillins




Cephalosporins



Cephamycins




Carbapenems




Monobactam

Gram-Negative Bacterial Cell



β -Lactam Activity Against Gram-Negatives




Penicillin-Binding
Proteins (PBPs)

Cartoon courtesy of L. Westblade


Periplasmic
space

Cell wall
(peptidoglycan)⁹


β -Lactam Activity Against Gram-Negatives



Penicillin-Binding Proteins (PBPs)




β -Lactam Resistance in Gram-negatives: β -Lactamases



Cartoon courtesy of L. Westblade

β -Lactam Resistance in Gram-negatives: Porin Channel Obstructed or Efflux Pump Altered




Outer membrane

Periplasmic space

Cell wall (peptidoglycan)²


Penicillin-Binding Proteins (PBPs)



β -Lactamases Produced by Gram-negatives

- ◆ Enzymes that hydrolyze the β -lactam ring, inactivating the β -lactam
- ◆ Hundreds of different types including:
 - Penicillinases
 - ESBLs
 - Carbapenemases
 - AmpCs
- ◆ Selectively inactivate various β -lactam antimicrobial agents


Mechanism of Action of β -Lactamases



**β -Lactam
(Carbapenem)**

**β -lactamase
(Carbapenemase)**





**Hydrolyzed
 β -Lactam
(Carbapenem)**

Concepts Related to AST

- ◆ When an enzyme produced by a bacterium **hydrolyzes** an antimicrobial agent, the enzyme inactivates that agent
- ◆ **Low level activity** of an enzyme implies that the MICs are reduced to a lesser extent than occurs with enzymes that exhibit higher level activity
 - MICs may be in the resistant, intermediate or even the high end of the susceptible range
 - Carbapenemases can exhibit low level or higher level activity

Carbapenem-Resistant Enterobacteriaceae

- ◆ Two mechanisms may result in carbapenem MICs or zone diameters in the “I” or “R” range among Enterobacteriaceae
 - Carbapenemase production
 - Cephalosporinase or ESBL together with porin loss
 - Some AmpC β -lactamases and ESBLs have low-level carbapenem-hydrolyzing activity
 - Porin loss limits entry of the carbapenem into the cell

GNR β -Lactamases (Non-carbapenemase)

Class	Examples	Produced by:	Notes
A	ESBLs [TEM, SHV, CTX-M]	<i>K. pneumoniae</i> and other Enterobacteriaceae	Most inhibited by β -lactamase inhibitors Usually plasmid-mediated; Can confer carbapenem resistance if other “R” mechanisms are present (e.g., porin modification)
B	---	---	---
C	AmpC	Enterobacteriaceae and some non-fermenters	Inducible in some genera (SPACE/SPICE organisms); Not inhibited by clavulanic acid
D	---	---	---

Adapted from Queenan & Bush. 2007. Clin Microbiol Rev. 20:440.

Bush & Jacoby. 2010. AAC. 54:969; Bush, K. 2013. Ann NY Acad Sci 1277:84.

GNR Carbapenemases

Class	Examples	Produced by:	Notes
A	Serine carbapenemases: KPC, SME	<i>K. pneumoniae</i> and other Enterobacteriaceae <i>S. marcescens</i>	Usually plasmid- mediated (not SME)
B	MBL carbapenemases: e.g. NDM, VIM, IMP, GIM, SPM	<i>P. aeruginosa</i> Enterobacteriaceae <i>Acinetobacter</i> <i>S. maltophilia</i>	Inhibited by EDTA Do not hydrolyze aztreonam
C	---	---	---
D	OXA carbapenemases	<i>Acinetobacter baumannii</i> Pseudomonads Enterobacteriaceae	Weakly hydrolyze carbapenems

Adapted from Queenan & Bush. 2007. Clin Microbiol Rev. 20:440.

Bush & Jacoby. 2010. AAC. 54:969; Bush, K. 2013. Ann NY Acad Sci 1277:84.

**Most Common
Carbapenemases**

KPC

(*Klebsiella pneumoniae* Carbapenemase)

- ◆ Most common carbapenemase in USA
- ◆ First report 1996 from North Carolina
- ◆ Usually a high level of enzyme can be produced
- ◆ Mostly *K. pneumoniae*, also *K. oxytoca*, *E. coli*, *C. freundii*, *Enterobacter* spp., *Salmonella*, *Serratia* spp., *P. aeruginosa* and other GNRs
- ◆ Plasmid with KPC gene generally has other R genes including genes for ESBLs


Metallo β -Lactamase (MBL) Carbapenemase

- ◆ **NDM** (New Delhi MBL) is the most common MBL worldwide; frequently encountered in **India and Pakistan**
- ◆ First report **2008** in a Swedish patient who was hospitalized in **India**
- ◆ Called MBL because **zinc** is required for activity
- ◆ Mostly *K. pneumoniae* and *E. coli*
- ◆ *bla*_{NDM} gene is highly mobile
- ◆ Also includes *bla*_{IMP}, *bla*_{VIM}

OXA Carbapenemase

- ◆ First described in *Acinetobacter baumannii* in 1985
- ◆ OXA-48
 - Commonly found in Europe and Africa; relatively rare in USA
 - First reported in 2008 in Turkey
 - Mostly *K. pneumoniae*, *E. coli*
- ◆ Many OXA-48-like variants described to date (OXA-181, OXA-232)
- ◆ Weakly hydrolyze carbapenems and cephalosporins (OXA-48 has greater hydrolytic properties than some other OXAs against carbapenems)

Number of MBL-producing Gram-Negative Rods Reported in the U.S. to the CDC (as of 4/27/2016)



Slide courtesy of
B. Limbago, CDC

Carbapenemase Spread and Diversity

- ◆ A single patient may harbor several species containing *bla*_{KPC}
 - One patient had KPC-producing *K. pneumoniae*, *E. coli*, and *Serratia marcescens*
Sidjabat et al. 2009. Clin Infect Dis. 49:1736.
 - ◆ A single patient may harbor several species with different carbapenemases
 - One UCLA patient
 - KPC (*K. pneumoniae*)
 - SME (*S. marcescens*)
- Pollett et al. 2014. J Clin Microbiol. 52:4003.

Breakpoints

Breakpoint Reminders

- ◆ CLSI and FDA set / revise breakpoints in USA
- ◆ Carbapenem breakpoints for Enterobacteriaceae were updated by:
 - CLSI in 2010
 - FDA in 2012 (ertapenem, imipenem); 2013 (meropenem)
 - Note: for doripenem, FDA lists “S” only breakpoint of ≤ 0.5 $\mu\text{g/ml}$ (zone ≥ 23 mm); 2012
- ◆ Current FDA breakpoints may not be updated on commercial systems
 - Clinical laboratory must perform a **verification** for a commercial system if using breakpoints other than those that are FDA-cleared on that system

Check with manufacturer of your commercial system for breakpoint status.

Where are the current CLSI breakpoints?

M100S 26th ed Table 2A (January 2016)

Table 2A-1. Zone Diameter and Minimal Inhibitory Concentration Interpretive Standards for *Enterobacteriaceae*

Testing Conditions		Routine QC Recommendations (See Tables 4A and 5A for acceptable QC ranges.)	
Medium:	Disk diffusion: MHA Broth dilution: CAMHB Agar dilution: MHA	<i>Escherichia coli</i> ATCC® 25922	
Inoculum:	Growth method or direct colony suspension, equivalent to a 0.5 McFarland standard	<i>Pseudomonas aeruginosa</i> ATCC® 27853 (for carbapenems)	
Incubation:	35°C ± 2°C; ambient air Disk diffusion: 16–18 hours Dilution methods: 16–20 hours	<i>Escherichia coli</i> ATCC® 35218 (for β-lactam/β-lactamase inhibitor combinations)	
		When a commercial test system is used for susceptibility testing, refer to the manufacturer's instructions for QC test recommendations and QC ranges.	

* ATCC® is a registered trademark of the American Type Culture Collection.

Refer to Tables 3A, 3B, and 3C for additional testing recommendations, reporting suggestions, and QC.

General Comments

- (1) For disk diffusion, test a maximum of 12 disks on a 150-mm plate and no more than 6 disks on a 100-mm plate; disks should be placed no less than 24 mm apart, center to center (see M02-A12, Subchapter 3.6). Each zone diameter should be clearly measurable; overlapping zones prevent accurate measurement. Measure the diameter of the zones of complete inhibition (as judged by the unaided eye), including the diameter of the disk. Hold the Petri plate a few inches above a black background illuminated with reflected light. The zone margin should be considered the area showing no obvious, visible growth that can be detected with the unaided eye. Ignore faint growth of tiny colonies that can be detected only with a magnifying lens at the edge of the zone of inhibited growth. Strains of *Proteus* spp. may swarm into areas of inhibited growth around certain antimicrobial agents. With *Proteus* spp., ignore the thin veil of swarming growth in an otherwise obvious zone of growth inhibition. With trimethoprim and the sulfonamides, antagonists in the medium may allow some slight growth; therefore, disregard slight growth (20% or less of the lawn of growth) and measure the more obvious margin to determine the zone diameter.
- (2) When fecal isolates of *Salmonella* and *Shigella* spp. are tested, only ampicillin, a fluoroquinolone, and trimethoprim-sulfamethoxazole should be reported routinely. In addition, for extraintestinal isolates of *Salmonella* spp., a third-generation cephalosporin should be tested and reported, and chloramphenicol may be tested and reported if requested. Susceptibility testing is indicated for typhoidal *Salmonella* (*S. Typhi* and *Salmonella* Paratyphi A–C) isolated from extraintestinal and intestinal sources. Routine susceptibility testing is not indicated for nontyphoidal *Salmonella* spp. isolated from intestinal sources. **In contrast, susceptibility testing is indicated for all *Shigella* isolates.**
- (3) The dosage regimens shown in the comments column below are those needed to achieve plasma drug exposures (in adults with normal renal and hepatic functions) on which breakpoints were based. When implementing new breakpoints, it is strongly recommended that laboratories share this information with infectious diseases practitioners, pharmacists, pharmacy and therapeutics committees, and infection control committees.



CLINICAL AND
LABORATORY
STANDARDS
INSTITUTE®

26th Edition

M100S

Performance Standards for Antimicrobial Susceptibility Testing

➤

This document provides updated tables for the Clinical and Laboratory Standards Institute antimicrobial susceptibility testing standards M02-A12, M07-A10, and M11-A8.

An informational supplement for global application developed through the Clinical and Laboratory Standards Institute consensus process.

Subcommittee on Antimicrobial Susceptibility Testing
This document is protected by international copyright laws.

27

Where are the current FDA breakpoints?

U.S. Department of Health and Human Services

U.S. Food and Drug Administration
Protecting and Promoting Your Health

A to Z Index | Follow FDA | En Español

Search FDA

Home | Food | Drugs | Medical Devices | Radiation-Emitting Products | Vaccines, Blood & Biologics | Animal & Veterinary | Cosmetics | Tobacco Products

About FDA

Home > About FDA > FDA Organization > Office of Medical Products and Tobacco > About the Center for Drug Evaluation and Research

About the Center for Drug Evaluation and Research

CDER Offices and Divisions

Drug Safety Oversight Board

Jobs at the Center for Drug Evaluation and Research (CDER)

Meeting Presentations (Drugs)

CDER Exclusivity Board

FAQs about CDER

Reports & Budgets (CDER)

Manual of Policies & Procedures (CDER)

Contact CDER

Antibacterial and Antifungal Product Labeling: Microbiology Susceptibility Interpretive Criteria (Breakpoints) and Quality Control Parameter Updates

SHARE | TWEET | LINKEDIN | PIN IT | EMAIL | PRINT

Table Last Updated: January 21, 2016
New or Changes are in bold.

This table lists approved Antibacterial and Antifungal products whose Microbiology subsection, Susceptibility Interpretive Criteria and Quality Control Parameter information in the labeling, has been reviewed and is current. The labels for these products are posted on the [Drugs@FDA](#) or [Daily Med](#) websites.

Active Ingredient	Product Name	Company Name	Type of Application	Application Number	Date of Most Recent FDA Review of Microbiology Susceptibility Interpretive Criteria*
Amoxicillin	Amoxil (amoxicillin) Chewable Tablets	Dr. Reddy's Laboratories	NDA	50-542	09/24/15

Updated
1/21/16

<http://www.fda.gov/AboutFDA/CentersOffices/OfficeofMedicalProductsandTobacco/CDER/ucm275763.htm>

Google: "FDA Interpretive Criteria Updates" and you will see list of dates of updated breakpoints

To view the FDA breakpoints for a specific antimicrobial agent, click on: “Dailymed” or “Drugs@FDA”

U.S. Department of Health and Human Services

FDA U.S. Food and Drug Administration
Protecting and Promoting *Your Health*

A to Z Index | Follow FDA | En Español

Search FDA

Home | Food | Drugs | Medical Devices | Radiation-Emitting Products | Vaccines, Blood & Biologics | Animal & Veterinary | Cosmetics | Tobacco Products

About FDA

Home > About FDA > FDA Organization > Office of Medical Products and Tobacco > About the Center for Drug Evaluation and Research

About the Center for Drug Evaluation and Research

- CDER Offices and Divisions
- Drug Safety Oversight Board
- Jobs at the Center for Drug Evaluation and Research (CDER)
- Meeting Presentations (Drugs)
- CDER Exclusivity Board
- FAQs about CDER
- Reports & Budgets (CDER)
- Manual of Policies & Procedures (CDER)
- Contact CDER

Antibacterial and Antifungal Product Labeling: Microbiology Susceptibility Interpretive Criteria (Breakpoints) and Quality Control Parameter Updates

SHARE | TWEET | LINKEDIN | PIN IT | EMAIL | PRINT

Table Last Updated: January 21, 2016
New or Changes are in bold.

This table lists approved Antibacterial and Antifungal products whose Microbiology subsection, Susceptibility Interpretive Criteria and Quality Control Parameter information in the labeling, has been reviewed and is current. The labels for these products are posted on the [Drugs@FDA](#) or [Daily Med](#) websites.

Active Ingredient	Product Name	Company Name	Type of Application	Application Number	Date of Most Recent FDA Review of Microbiology Susceptibility Interpretive Criteria*
Amoxicillin	Amoxil (amoxicillin) Chewable Tablets	Dr. Reddy's Laboratories	NDA	50-542	09/24/15

Type in antimicrobial agent name, then click on “Label Information” for the antimicrobial agent, which will open up a pdf of the package insert.

The image shows a screenshot of the DailyMed website, which is part of the U.S. National Library of Medicine. The page features a search bar with the placeholder text "Enter drug, NDC code, drug class, or Set ID". Below the search bar, there are navigation tabs for "ALL DRUGS", "HUMAN DRUGS", and "ANIMAL DRUGS". The "Label Information" tab is highlighted. The page also includes a "NEWS" section with a post dated March 7, 2016, and a "Drug Approval Reports by Month" section. The FDA logo and "U.S. Food and Drug Administration" text are visible in the top right corner of the screenshot.

U.S. NATIONAL LIBRARY OF MEDICINE

REPORT ADVERSE EVENTS | RECALLS

DAILYMED

ALL DRUGS | HUMAN DRUGS | ANIMAL DRUGS

Enter drug, NDC code, drug class, or Set ID

MORE WAYS TO SEARCH: [ADVANCED SEARCH](#) [BROWSE DRUG CLASSES](#) [LABEL INFORMATION](#)

This website contains 84522 drug listings as submitted to the Food and Drug Administration. At the present time, this Web site does not contain a complete listing of labels for approved drugs.

NEWS

[DailyMed Announcements](#)

Posted: March 7, 2016
DailyMed Now Served Over HTTPS Only
DailyMed has been updated to serve all content over HTTPS only.

The base URL for the web application is: <https://dailymed.nlm.nih.gov/>

The base URL for the web service is:
<https://dailymed.nlm.nih.gov/dailymed/services/>

FDA U.S. Food and Drug Administration
Protecting and Promoting Your Health

A to Z Index | Follow FDA | En Español

SEARCH

Most Popular Searches

Home Food Drugs Medical Devices Radiation-Emitting Products Vaccines, Blood & Biologics Animal & Veterinary Cosmetics Tobacco Products

FDA Home Drug Databases Drugs@FDA

Drugs@FDA
FDA Approved Drug Products

FAQ | Instructions | Glossary | Contact Us
Drugs@FDA Demo | What's New in Drugs@FDA

Search by Drug Name, Active Ingredient, or Application Number

Enter at least three characters: [Advanced Search](#)

Browse by Drug Name

A B C D E F G H I J K L M N O P Q R S T U V W X Y Z 0-9

Drug Approval Reports by Month

Disclaimer

FDA/Center for Drug Evaluation and Research
Office of Communications
Division of Online Communications
Update Frequency: Daily

Then Go to “Clinical Pharmacology section” in the package insert (subsection Microbiology) for FDA BPs

- 11 DESCRIPTION**
- 12 CLINICAL PHARMACOLOGY**
 - 12.1 Mechanism of Action
 - 12.2 Pharmacodynamics
 - 12.3 Pharmacokinetics
 - 12.4 Microbiology
- 13 NON-CLINICAL TOXICOLOGY**
 - 13.1 Carcinogenesis, Mutagenesis, and Impairment of Fertility

Table 6: Susceptibility Test Result Interpretive Criteria for Doripenem

Pathogen	Minimum Inhibitory Concentrations (µg/mL)	Disk Diffusion (zone diameters in mm)
	Susceptible*	Susceptible*
<i>Enterobacteriaceae</i>	≤ 0.5	≥ 23
<i>Pseudomonas aeruginosa</i>	≤ 2	≥ 24
<i>Acinetobacter baumannii</i>	≤ 1	≥ 17
<i>Streptococcus anginosus</i> group (<i>S. constellatus</i> and <i>S. intermedius</i>)	≤ 0.12	≥ 24
Anaerobes	≤ 1	n/a

* The current absence of resistant isolates precludes defining any results other than "Susceptible". Isolates yielding MIC or disk diffusion results suggestive of "Nonsusceptible" should be subjected to additional testing.
n/a = not applicable

FDA Breakpoints for Doripenem

Where can I find information about “verification” of ASTs?

ASM Press

Browse by Category

ASMPRESS Categories

- Applied and Industrial Microbiology
- ASM Meetings - Abstract CDRoms
- Bacterial Pathogenesis
- Catalog
- Clinical Microbiology
- Cumitechs - Practical Clinical Guidelines
- E-BOOK INFORMATION
- eBooks iPad App
- Environmental Microbiology
- Food Microbiology
- Fungi and Fungal Pathogenesis
- General Interest
- History of Science
- Immunology
- Medical Microbiology, Pathogenesis
- Microbial Genetics and Molecular Biology
- Molecular Biology, Genetics
- New and Forthcoming Titles
- Online and Apps
- Online Periodicals

Pricing

List and ASM member price (WEB DOWNLOAD) \$24.95 **Buy**

Cumitech 31A: Verification and Validation of Procedures in the Clinical Microbiology Laboratory (WEB DOWNLOAD)

Access Period: 14 days
Author: Coordinating Editor: Susan E. Sharp
Book ISBN or Item Number: 978-1-55581-530-1

RightsLink Permissions

Coordinating Editor
Susan E. Sharp

Clinical Microbiology
NEWSLETTER

CMN Stay Stay

CMN
Vol. 35, No. 13
July 1, 2013
www.cmnewsletter.com

Verification of Antimicrobial Susceptibility Testing Methods: a Practical Approach

Jean B. Patta,¹ Susan Sharp,² and Susan Novak-Weekley,³ ¹Division of Healthcare Quality Protection for Disease Control and Prevention, Atlanta, Georgia, ²Northwest Permanente Physicians and Staff Oregon, and ³SCPMG Regional Reference Laboratories, North Hollywood, California

IN THIS ISSUE

103 Verification of Antimicrobial Susceptibility Testing Methods: a Practical Approach

109 Case Report: Fatal Pulmonary *Mycobacterium abscessus* Infection in an Immunosuppressed Patient

Abstract

The process of verifying an antimicrobial susceptibility testing (AST) system can be very complex. There are many different reasons why verification might be necessary, such as implementing a new method in the laboratory or implementing non-FDA interpretive criteria or breakpoints on an FDA-cleared AST system. The Clinical Laboratory Improvement Amendment (CLIA) provides some general guidance, but ultimately, it is the responsibility of a laboratory director to decide the composition of a verification study protocol. Variables to consider are what methods should be compared, what and how many isolates should be tested, how the results will be compared, and what study results will result in an acceptable study outcome. This article provides some general guidelines for developing and conducting a verification study of an AST system.

Introduction

The Clinical Laboratory Improvement Amendment (CLIA) regulation requires laboratories to verify the performance of their testing methods. When laboratories implement revised breakpoints face a significant obstacle, because only FDA breakpoints can be used on FDA-approved devices. Use of alternative breakpoints requires an in-house verification

CLINICAL AND LABORATORY STANDARDS INSTITUTE

1st Edition

M52

Verification of Commercial Microbial Identification and Antimicrobial Susceptibility Testing Systems


This guideline includes recommendations for verification of commercial US Food and Drug Administration–cleared microbial identification and antimicrobial susceptibility testing systems by clinical laboratory professionals to fulfill regulatory or quality assurance requirements for the use of these systems for diagnostic testing.

A guideline for US application developed through the Clinical and Laboratory Standards Institute consensus process.

Revised by: CLSIA Medical Center
 This document is protected by copyright. CLSI order # 0. Downloaded on 08/27/2015

AST Verification Requirements

- ◆ In USA, laboratories must perform verification study **when implementing a new test or new breakpoints**
 - CLIA requirement
- ◆ Must perform if using:
 - Disk diffusion
 - Commercial system that is FDA-cleared for the new breakpoints
 - Commercial system that is not FDA-cleared for the new breakpoints



FDA-cleared AST
...Prior to patient testing

AST In-House Verification Study - Example*

Step	Details
Obtain Isolates N ≈ 30	<ul style="list-style-type: none">• Obtain isolates from a reference laboratory or other source (“reference set”); “Reference set”:<ul style="list-style-type: none">- is provided with meropenem MIC and S, I, R results that were obtained from method verified with current CLSI/FDA BPs- includes both “S” and “R” phenotypes with variety of MICs
Test by Automated MIC	<ul style="list-style-type: none">• Evaluate carbapenems routinely tested in your laboratory• Interpret MICs (M100S 26th ed)
Review MIC and S, I, R results	<ul style="list-style-type: none">• Compare MIC results from automated AST system to “reference” MIC results (essential agreement)• Compare S, I, R results from automated AST system to “reference” S, I, R results (categorical agreement)• Some results may not agree; need plan to arbitrate (e.g., send to another lab that uses CLSI reference MIC method)
Write verification report	<ul style="list-style-type: none">• Submit for laboratory director approval

* **verify current CLSI/FDA meropenem BPs on commercial automated AST system using set of organisms obtained from a reference laboratory**

Example: Meropenem In House Verification Study

Test Method: Automated AST System

Reference (Comparator): “Reference Set” from Lab B

Organism	Isolates No. (No. R)	EA		CA		VME		ME	
		#	%	#	%	#	%	#	%
<i>E. coli</i> 2 CR (KPC)	10 (2)	10	100	10	100	0	0	0	0
<i>K. pneumoniae</i> 5 CR (4 KPC, 1 NDM)	10 (5)	9	90	9	90	1	20*	0	0
Other Enterobacteriaceae (1 SME, 1 NDM)	10 (2)	9	90	10	100	0	0	0	0
Total	30 (9)	28	93.3	29	97	1	11*	0	0

CR = carbapenem resistant; VME = very major error; ME = major error

* Unacceptable result → repeat, send to second reference lab for additional testing.
Isolates may have lost plasmid. Avoid excessive subbing to avoid plasmid loss.

Enterobacteriaceae

Interpreting Carbapenem Results

**The best way to detect CRE is to use
current breakpoints!**

CLSI M100S 26th ed Table 2A

Enterobacteriaceae and Carbapenems

regimen of 1 g every 8 h.
See comment (7).

CARBAPENEMS

(25) Following evaluation of PK-PD properties, limited clinical data, and MIC distributions that include recently described carbapenemase-producing strains, revised interpretive criteria for carbapenems were first published in June 2010 (M100-S20-U) and are listed below. Because of limited treatment options for infections caused by organisms with carbapenem MICs or zone diameters in the intermediate range, clinicians may wish to design carbapenem dosage regimens that use maximum recommended doses and possibly prolonged intravenous infusion regimens, as has been reported in the literature.¹⁻⁴ Consultation with an infectious diseases practitioner is recommended for isolates for which the carbapenem MICs or zone diameter results from disk diffusion testing are in the intermediate or resistant ranges.

Laboratories using *Enterobacteriaceae* MIC interpretive criteria for carbapenems described in M100-S20 (January 2010) should perform the MHT, the Carba NP test, and/or a molecular assay when isolates of *Enterobacteriaceae* are suspicious for carbapenemase production based on imipenem or meropenem MICs of 2–4 µg/mL or ertapenem MIC of 2 µg/mL (refer to Tables 3B and 3C). After implementation of the current interpretive criteria, the MHT does not need to be performed other than for epidemiological or infection control purposes (refer to Table 3B).

The following information is provided as background on carbapenemases in *Enterobacteriaceae* that are largely responsible for MICs and zone diameters in the intermediate and resistant ranges, and thus the rationale for setting revised carbapenem breakpoints:

- The clinical effectiveness of carbapenem treatment of infections produced by isolates for which the carbapenem MIC or disk diffusion test results are within the intermediate range is uncertain due to lack of controlled clinical studies.
- Imipenem MICs for *Proteus* spp., *Providencia* spp., and *Morganella morganii* tend to be higher (eg, MICs in the intermediate or resistant range) than meropenem or doripenem MICs. These isolates may have elevated imipenem MICs by mechanisms other than production of carbapenemases.

B	Doripenem	10 µg	≥23	–	20–22	≤19	≤1	–	2	≥4	(26) Interpretive criteria are based on a dosage regimen of 500 mg every 8 h.
B	Ertapenem	10 µg	≥22	–	19–21	≤18	≤0.5	–	1	≥2	(27) Interpretive criteria are based on a dosage regimen of 1 g every 24 h.
B	Imipenem	10 µg	≥23	–	20–22	≤19	≤1	–	2	≥4	(28) Interpretive criteria are based on a dosage regimen of 500 mg every 6 h or 1 g every 8 h.
B	Meropenem	10 µg	≥23	–	20–22	≤19	≤1	–	2	≥4	(29) Interpretive criteria are based on a dosage regimen of 1 g every 8 h.

CLSI Carbapenem Breakpoints / Dosage Comments

S

I

R

S

I

R

B	Doripenem	10 µg	≥23	-	20-22	≤19	≤1	-	2	≥4	(26) Interpretive criteria are based on a dosage regimen of 500 mg every 8 h.
B	Ertapenem	10 µg	≥22	-	19-21	≤18	≤0.5	-	1	≥2	(27) Interpretive criteria are based on a dosage regimen of 1 g every 24 h.
B	Imipenem	10 µg	≥23	-	20-22	≤19	≤1	-	2	≥4	(28) Interpretive criteria are based on a dosage regimen of 500 mg every 6 h or 1 g every 8 h.
B	Meropenem	10 µg	≥23	-	20-22	≤19	≤1	-	2	≥4	(29) Interpretive criteria are based on a dosage regimen of 1 g every 8 h.

(26) Interpretive criteria are based on a dosage regimen of 500 mg every 8 h.

(27) Interpretive criteria are based on a dosage regimen of 1 g every 24 h.

(28) Interpretive criteria are based on a dosage regimen of 500 mg every 6 h or 1 g every 8 h.

(29) Interpretive criteria are based on a dosage regimen of 1 g every 8 h.

CLSI Carbapenem Dosage Comment

“Because of limited treatment options for infections caused by organisms with carbapenem MICs or zone diameters in the intermediate range, clinicians may wish to design carbapenem dosage regimens that use maximum recommended doses and possibly prolonged intravenous infusion regimens as has been reported in the literature.”

Maintain higher antimicrobial agent levels in patients for longer periods of time.

Enterobacteriaceae - Carbapenem Breakpoints (MIC $\mu\text{g/ml}$)¹

Agent	Current CLSI & FDA Breakpoints			Old CLSI Breakpoints		
	Susc	Int	Res	Susc	Int	Res
Ertapenem	≤ 0.5	1	≥ 2	≤ 2	4	≥ 8
Imipenem	≤ 1	2	≥ 4	≤ 4	8	≥ 16
Meropenem	≤ 1	2	≥ 4	≤ 4	8	≥ 16

Agent	Current CLSI Breakpoints			Old CLSI Breakpoints
Doripenem ²	≤ 1	2	≥ 4	none

¹CLSI M100 26th ed; also lists corresponding disk diffusion breakpoints

²FDA breakpoint "S" only = ≤ 0.5

Carbapenemase

Testing for the Enzymes

Do we need to determine if a carbapenemase is present in CRE?

◆ Patient care

- No, when using the “current” CLSI breakpoints
- Yes, if using “old” CLSI breakpoints

◆ Infection control

- Yes, if outbreak suspected; possibly in other settings

◆ Epidemiology / research

- Yes, to better understand emerging resistance and plan for “challenges”

Some say MUST identify resistance mechanism in all CRE as carbapenemase-producers are more worrisome than non-carbapenemase-producing CRE.

Introduction to Tables 3B and 3C. Tests for Carbapenemases in Enterobacteriaceae, *Pseudomonas aeruginosa*, and *Acinetobacter* spp.


	MHT	Carba NP	Molecular
Use	Enterobacteriaceae	Enterobacteriaceae <i>P. aeruginosa</i> <i>Acinetobacter</i>	Enterobacteriaceae <i>P. aeruginosa</i> <i>Acinetobacter</i>
Strengths	Simple	Rapid	Determines type of carbapenemase
Limitations	Some false pos (eg, ESBL/AmpC + porin) Some false neg (eg NDM) Enterobacteriaceae only	Special “fresh” reagents Some invalid results False neg for OXA-type carbapenemase	Special reagents Specific to targeted gene Difficult for some micro labs to implement

Modified Hodge Test for Carbapenemases

Table 3B. The Modified Hodge Confirmatory Test for Suspected Carbapenemase Production in *Enterobacteriaceae*

NOTE: If using FORMER minimal inhibitory concentration (MIC) interpretive criteria for carbapenems described in M100-S20 (January 2010), please refer to modifications in Table 3B-1 below.


Test	Confirmatory Test												
When to Do This Test:	For epidemiological or infection control purposes. NOTE: No change in the interpretation of carbapenem susceptibility test results is required for carbapenemase-positive isolates.												
Test Method	MHT												
Medium	MHA												
Antimicrobial Concentration	Ertapenem disk 10 µg or Meropenem disk 10 µg												
Inoculum	<p>(1) Prepare a 0.5 McFarland standard suspension (using either direct colony suspension or growth method) of <i>E. coli</i> ATCC® 25922 (the indicator organism) in broth or saline, and dilute 1:10 in saline or broth. Inoculate an MHA plate as for the routine disk diffusion procedure. Allow the plate to dry 3 to 10 minutes. Place the appropriate number of ertapenem or meropenem disks on the plate as noted below and shown in Figures 1 and 2.</p> <p>(2) Using a 10-µL loop or swab, pick 3 to 5 colonies of test or QC organism grown overnight in broth or saline, and streak a straight line out from the edge of the disk. The streak should be at least 20–25 mm in length and extend to the center of the disk as noted below and shown in Figures 1 and 2.</p> <p>Capacity of small and large MHA plates (100-mm or 150-mm diameter, respectively):</p> <table border="1"> <thead> <tr> <th></th> <th>Small</th> <th>Large</th> </tr> </thead> <tbody> <tr> <td>Disks</td> <td>1</td> <td>1–4</td> </tr> <tr> <td>Test isolates</td> <td>1</td> <td>1–6</td> </tr> <tr> <td>QC isolates</td> <td>2</td> <td>2</td> </tr> </tbody> </table>		Small	Large	Disks	1	1–4	Test isolates	1	1–6	QC isolates	2	2
	Small	Large											
Disks	1	1–4											
Test isolates	1	1–6											
QC isolates	2	2											
Incubation Conditions	35°C ± 2°C; ambient air												
Incubation Length	16–20 hours												



M100S 26th ed. Table 3B.

Carba NP Test for Carbapenemase Production

- ◆ Isolated colonies (lyse)
- ◆ Hydrolysis of imipenem
- ◆ Detected by change in pH of indicator (red changes to yellow/orange)
- ◆ Rapid <2h
- ◆ Microtube method




NO +
imipenem imipenem

Nordmann et al. 2012. *Emerg Infect Dis.* 18:1503.
Tijet et al. 2013. *Antimicrob Agents Chemother.* 57:4578.
Vasoo et al. 2013. *J Clin Microbiol.* 51:3092.
Dortet et al. 2014. *J Med Microbiol.* 63:772.
Dortet et al. 2014. *Antimicrob Agents Chemother.* 58:2441.



Carba NP Test - Examples



Blank

Neg

KPC

OXA48

OXA181

NDM

IMP

VIM

SME

Pos

Invalid

Neg

Pos


Pos

Pos

Pos

Courtesy of Shaun Yang, P. Hemarajata


Modified Hodge Test - Examples



Neg Control



**Pos Control
KPC**



**NDM
False neg**




**OXA 232
Positive**



**SME
Positive**

Modified Hodge Test False Positive



CIM Test – for Carbapenemase (Carbapenem Inactivation Method)

◆ Principle:

- A suspension of bacteria is incubated with a standard meropenem disk.
- If the organism produces carbapenemase, it will inactivate the meropenem in the disk.
- Following 2 h incubation, the meropenem disk is removed from the suspension and placed on a lawn of *E. coli* ATCC 25922.
- Following overnight incubation, carbapenemase activity is demonstrated by a loss of meropenem activity in the disk.



#1 carbapenemase positive
#2 carbapenemase negative

Neo-Rapid CARB Screen Kit

- ◆ Commercial kit; similar to Carba NP
- ◆ Enterobacteriaceae and *P. aeruginosa*
- ◆ Tablets
 - Imipenem + indicator
 - Negative control
- ◆ ≤2 hours
- ◆ Research use only



Molecular Tests for Carbapenemases¹

- ◆ **Biofire**²
 - KPC
- ◆ **Nanosphere**²
 - KPC, NDM, OXA, IMP, VIM
- ◆ **Cepheid**³
 - KPC, NDM, OXA-48, IMP-1, VIM
- ◆ **BD Max**⁴
 - KPC, NDM, OXA-48
- ◆ **NucliSENS EasyQ KPC**⁴
 - KPC
- ◆ **Check-Points CPE**⁴
 - KPC, NDM, OXA-48, OXA-181, IMP, VIM

¹ Not all inclusive

² FDA-cleared; for use with positive blood culture broth

³ FDA-cleared; for use with isolated colonies

⁴ RUO (Research Use Only)

CRE

Reporting Results

Testing and Reporting Additional Agents

Specimen: Blood
Diagnosis: Pneumonia
Klebsiella pneumoniae

Final Report with
Optional Comment

MIC ($\mu\text{g/ml}$)

amikacin	>32 R
cefepime	>32 R
ceftriaxone	>32 R
ciprofloxacin	>2 R
ertapenem	>4 R
gentamicin	>10 R
meropenem	4 R
piperacillin-tazo	>128/4 R
tobramycin	>10 R
trimeth-sulfa	>4/76 R

“This *Klebsiella pneumoniae* has unusual carbapenem results and is considered carbapenem-resistant Enterobacteriaceae (CRE); Infectious Diseases Consult suggested”

Treatment Options for CRE

- ◆ Polymyxins
 - Colistin, polymyxin B
- ◆ Tigecycline
- ◆ Minocycline
- ◆ Aminoglycosides (not tobramycin)
- ◆ Fosfomycin
- ◆ Ceftazidime-avibactam (class A producers only)

Usually a combination of antimicrobial agents (often with a polymyxin) are used

Supplemental Antimicrobial Agents to Consider for AST of CRE

Antimicrobial Agent	Enterobacteriaceae (CRE)
Minocycline	Yes
Tigecycline ¹	Yes (excluding <i>Pro/Prov/Morg</i>)
Colistin (or polymyxin B)	Yes ²
Fosfomycin	Yes ³
Ceftazidime-avibactam	Yes

¹ Not on urine isolates. Breakpoints ($\mu\text{g/ml}$): No CLSI; FDA ≤ 2 S, 4 I, ≥ 8 R; EUCAST ≤ 1 S, > 2 R

² Breakpoints ($\mu\text{g/ml}$): No CLSI; EUCAST ≤ 2 S, > 2 R

³ Breakpoints ($\mu\text{g/ml}$): CLSI ≤ 64 S, ≥ 256 R for *E. coli* only; EUCAST ≤ 32 S, > 32 R for all Enterobacteriaceae

Availability of ASTs for Supplemental Drugs for CRE

Antimicrobial	Automated Systems	Disk Diffusion	Etest	CLSI Reference MIC method ¹
Minocycline	Some	Yes	Yes	Yes
Tigecycline	Yes	Yes	Yes	Yes
Colistin (or polymyxin B) ²	No	Yes, but poor performance, RUO	Poor performance, RUO	Yes
Fosfomycin ^{3,4}	No	Yes	Yes	Yes
Ceftaz-avibactam	Some	Yes, RUO	Yes, RUO	Yes

¹ performed by few laboratories

² no FDA-cleared commercial methods

³ only test reliably by agar methods

⁴ FDA-cleared for *E. coli* and *E. faecalis* only

Thank you for your attention!

**Please check this website again for
additional materials on CRE and other
topics in AST!**

**Brought to you by the Outreach
Working Group of the CLSI
Subcommittee on Antimicrobial
Susceptibility Testing**