


27 March 2019

To: Recipients of M02, 13th ed. and/or M07, 11th ed.
From: Jennifer K. Adams, MT(ASCP), MSHA
Vice President, Standards and Quality
Subject: Correction

This notice is intended to inform users of a correction made to CLSI documents M02, *Performance Standards for Antimicrobial Disk Susceptibility Tests*, 13th ed. and M07, *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically*, 11th ed. The correction is described below and shown as colored and/or stricken text in the flow chart.

Appendix D. Quality Control Protocol Flow Charts:

In section D1, Quality Control Protocol: Conversion From Daily to Weekly Testing (20- or 30-Day Plan), a flow chart box incorrectly states “Applicable QC strains are tested for another 20 tests days (Subchapter 4.7.2.1).” The flow chart box has been corrected to read “Applicable QC strains are tested for another 10 test days (Subchapter 4.7.2.1).”



If you require any additional clarification regarding these corrections, please contact CLSI Customer Service (customerservice@clsi.org).

We appreciate your commitment to CLSI and regret any inconvenience.